

VACUUM REGULATORS

5 YEAR WARRANTY



CE
0543

**HASSLE
FREE
SUCTION**

3 YEAR WARRANTY



Digital Vacuum Regulator

Analog Vacuum Regulator

Amvex Vacuum Regulators are the most technically advanced in the health care industry. The regulators' design lets the facility change from analog to digital with ease; upgrade from continuous to intermittent with an intermittent module; and allows for easy serviceability and maintenance.

Amvex Digital Regulators are the first of their kind in the marketplace. They improve patient care by providing greater accuracy and an easy to read gauge as compared to conventional analog regulators.

Our impressive line of regulators include: Continuous/Intermittent (combo), Continuous 2 Mode and 3 Mode, Pediatric Continuous 2 Mode, Pediatric Continuous/Intermittent (combo), Continuous High Vac.



FEATURES AND BENEFITS

Analog

- Analog gauge designed for clarity and appearance
- Operating ranges: 0-160 mmHg, 0-300 mmHg, and 0-760 mmHg
- Upgrade to Intermittent or Digital in minutes
- Same housing for both Continuous and Intermittent regulator
- Strong and resilient ABS plastic casing
- Removable back plate and gasket for easy cleaning



VR-CIU2-F2A



VR-C3U2-F2D

Digital

Includes all of the analog features plus....

- Patented state of the art Digital Display Technology
- Easy to read display
- Digital readout with no moving parts provides longer product life
- Accuracy is within +/-1%

PART CONFIGURATOR

VR-XX YY - XX Y (Z)

Vacuum Regulator Type:	
Continuous/Intermittent - 0-300 Gauge:	CI
Pediatric Intermittent - 0-160 Gauge:	PI
Pediatric Continuous - 0-160 Gauge:	P2
Continuous - High 0-760 Gauge:	CH
Continuous 2 Mode (Off/Reg) 0-300 Gauge:	C2
Continuous 3 Mode (Off/Reg/Full) 0-300 Gauge:	C3

Color Code:	
USA:	U
ISO:	I

Display:	
Digital:	D
Analog:	A

Outlet Connection:	
Tubing Nipple:	T
DISS Male:	D
1/8" FNPT:	2

Optional Colors:	
Yellow:	Y
Purple:	P
Grey:	E
Blue:	B
Red:	R
Green:	G



Color Options

Inlet Connection:	
1/8" FNPT:	F2
Ohmeda Male:	OM
Chemtron Male:	CM
DISS Handtight:	DH
DISS Hex Nut:	DN
Puritan Male:	PM
Medstar Male:	MM
Oxequip Male:	XM
British Male:	BM
French Male:	FM
German Male:	GM
Australian Male:	AM
Japanese Male:	JM

Visit our website for the Vacuum Regulator Part Configurator



www.amvex.com

Toll Free: 1-866-GO-AMVEX (1-866-462-6839) • Tel: (905) 764-7736 • Fax: (905) 764-7743 • Website: www.amvex.com

DEC 04