

Extended Latch-Valve Assembly Ohmeda



FRONT



SIDE



General Specifications:

The Extended Latch-Valve Assembly shall accommodate a 3/4" to 1-1/4" wall thickness.

A large color coded Front plate shall be used for ease of gas identification and aesthetic appeal.

Extended Latch-Valve Assemblies are available for Ohmeda Compatible, DISS, Puritan-Bennet Compatible, Chemetron Compatible and Oxequip Compatible Outlets.

All Extended Latch-Valve Assemblies shall be cleaned and degreased for medical gas service, factory assembled and tested.


Flow rate exceeds the requirements of NFPA and CSA.

Extended Latch-Valve Assemblies are UL Listed and ETL Listed.

Amico products comply with NFPA-99.



Technical Specifications



Model Numbers

O-FASE-XX-L-GAS

The Letters "XX" Defines the Latch-Valve Assembly Type:

- QD = Ohmeda Compatible
- CH = Chemetron Compatible
- PB = Puritan-Bennett Compatible
- DI = DISS
- OX = Oxequip/Medstar Compatible

The Letter "L" Defines the Language:

- U = English (NFPA) OXY, AIR, VAC & WAG
- E = English (CSA/ISO)
- F = French (CSA/ISO)

The Word "GAS" defines the Type of Gas:

- | | |
|--|-------|
| Oxygen | = OXY |
| Nitrogen * | = NIT |
| Medical Air | = AIR |
| Carbon Dioxide ** | = CO2 |
| Nitrous Oxide | = N2O |
| Medvac | = VAC |
| Waste Anesthetic Gas Disposal (NFPA) | = WAG |
| Anesthetic Gas Scavenging System (CSA) | = AGS |

NOTE:

- * Nitrogen available in DISS
- ** Carbon Dioxide not available in Puritan-Bennett