


Simplex Manifold

Reserve Manifold



Inch
[mm]



Features:

- Rigid copper pigtails with check valve
- 3000 psi working pressure
- CGA adapter port for pressure switch installation
- Wall mounted bracket included
- Staggard type for maximizing wall space
- High pressure master isolation valve
- High flow capacity regulator
- Cleaned for oxygen service
- Simplex Manifold complies with NFPA-99

Model Numbers

M-SIMP-E-XX-GAS

Simplex reserve manifold.

Replace XX with the total amount of cylinders.
Replace GAS with the appropriate gas.

For example: 10 Cylinder Medical Air = M-SIMP-E-10-AIR

Cylinders	Overall Manifold Length	
2	26 [660]	
4	39 [991]	
6	52 [1321]	Medical Air
8	65 [1651]	Oxygen
10	78 [1981]	Carbon Dioxide
12	91 [2311]	Nitrous Oxide
14	104 [2642]	Nitrogen
16	117 [2972]	
18	130 [3302]	
20	143 [3632]	

M-PRSW-RES

Pressure Switch set @ 1,100 psi decreasing pressure.



Amico Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 71 East Industry Court, Deer Park, NY 11729, USA
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 Tel: 905.764.0800 | Fax: 905.764.0862