



"Your Compressed Gas Equipment Specialists"

C-Block Manifold For Level 3 Systems

SPECIFICATIONS:

The C-BLOCK Manifold designed by Certified Medical Sales is recommended for NFPA Level 3 applications.

The manifold includes:

- 2 "H" Cylinder 2 Stage Regulators
- 2 – 5' long Hose Assemblies
- Brass Manifold Block with Mounting Holes
- 75 Psig Pressure Relief Valve
- 2 One Way Flow Check Valves
- 1 Barksdale "Reserve in Use" Alarm Switch
- 0-100 Psig Pressure Gauge



MODEL NUMBERS

C-BLOCK-OXY – Oxygen Manifold with Dual Regulators

C-BLOCK-N20 – Nitrous Manifold with Dual Regulators

C-BLOCK-NIT – Nitrogen Manifold with Dual Regulators

