



"Your Compressed Gas Equipment Specialists"

Mighty-Max Manifold For Level 3 Systems

SPECIFICATIONS:

The C-BLOCK Manifold designed by Certified Medical Sales is recommended for NFPA Level 3 applications.

The manifold includes:

- 4 "H" Cylinder 2 Stage Regulators
- 4 – 5' long Hose Assemblies
- Nickel Finish Brass Manifold Block with Mounting Holes
- 75 Psig Pressure Relief Valve
- 4 One Way Flow Check Valves
- 1 Barksdale "Reserve in Use" Alarm Switch
- 0-100 Psig Pressure Gauge
- Reserve/Primary Regulator Paddles



