

AC400B

Electric Bed with Five Functions



Product Benefits

- Easily removable head and foot boards
- Electrical head elevation (max 75 °) and knee gatch (max 40 °)
- Max trend +12°, reverse trend -8°
- Fifth wheel for ease of maneuverability
- Instant CPR release lever
- 500 lb. weight capacity
- LED night light
- Bed bumpers for added protection



Technical Specifications

The Amico Care Hospital Bed Collection provides Caregivers the control they need in the recovery environment. The simplistic design and distinctive features bring efficiency to today's medical and surgical environment.

The Amico Care AC400B features an easily removable head and foot boards which provides Caregivers easy access in the event of an emergency. The head and foot boards feature an exclusive blend of polypropylene for uncomplicated cleanup and durability.

The functions of the bed are electrically controlled and equipped with a back up battery that allows Caregivers to make adjustments to the bed during transport and unexpected power outages.

Dimensions		
Patient Surface	36" X 76"	(91 cm X 193 cm)
Overall Dimensions	42" X 83"	(107 cm X 211 cm)
Bed Frame Clearance	6"	(15 cm)
Bed Height		
Low Position	20"	(50 cm)
High Position	30"	(76 cm)
Bed Angles		
Max Head Elevation	75 °	
Knee Gatch	40°	
Max Trend / Reverse Trend	+12° / -8°	
Weight		
Overall Net Weight of Bed	243 lbs.	(110 kg)
Max Weight Limit	500 lbs	(227 kg)
Electrical		
Battery back-up(High low cycles)	2 - 12 V DC	
Voltage(please specify)	110 V/60 Hz or 220 V /50 Hz	
Caster & Braking System		
Caster Size	5"	(13 cm)

