

# Bedside Cabinets



BC-JD3-0BB-1CCL



BC-JDD-0XB-1CHM

## Jasper Series:

- High pressure laminate top and drawer fronts
- Casters (exposed or recessed) or self-levelling glides
- Different drawer configurations
- Seamless spill guard top available

## Oliver Series:

- Thermofoil top & drawer fronts
- Seamless for infection control
- Casters (exposed or recessed) or self-levelling glides
- Different drawer configurations
- Seamless spill guard top available



BC-ODD-1BT-2GCP



BC-OD3-1BA-4CHM



BC-LDD-0B-C1MC-L



BC-LDO-0B-G1HM-L

## Lucas Series:

- Anodized aluminum legs for durability
- Laminate or thermofoil tops and drawer fronts
- Casters (exposed) or self-levelling glides
- Different drawer configurations
- Seamless spill guard top available

# Customize Your Bedside Cabinets...

## Drawer Configurations



3 Drawers



Drawer/Door



Drawer/Open

## Features and Options


### Features



Steel side drawer construction



Vented back



Casters  
(2 locking, 2 unlocking)



Seamless spill guard top



Field replaceable drawer and door fronts

### Drawer Fronts



\*Bevelled



\*Convex



\*Shaker



\*Parallel



Flat

\*Only available in Thermofoil Material

### Hardware



Modern stainless steel



Stainless steel knob



Ornate stainless steel



Modern brass



Brass knob

### Finishes




Hard Rock Maple (HM)



Candlelight (CL)



Medium Oak (MO)



Medium Cherry (MC)



Chocolate Pear (CP)

Custom colors are available, please contact us for details

## Dimensions

Jasper Series	Height: 32"	Width: 18.5"	Depth: 18.75"
Oliver Series	Height: 32"	Width: 18.5"	Depth: 18.75"
Lucas Series	Height: 32"	Width: 19.5"	Depth: 19.5"