

AC200B

Electric Bed with Three Functions



Product Benefits

- Removable head and foot boards
- Instant CPR release lever
- Electrical head elevation (max 75 °) and knee gatch (max 40 °)
- 500 lbs. weight capacity



Technical Specifications

The AC200B bed provides Caregivers the control they need in the recovery environment. The simplistic design and distinctive features bring efficiency to today's medical and surgical environment.

The Amico Care AC200B features an open architecture for surfaces and removable head and foot boards which provides Caregivers easy access in the event of an emergency. The head and foot boards feature an exclusive blend of polypropylene for cleanup and durability.

The functions of the bed are electrically controlled and equipped with a back up battery that allows Caregivers to make adjustments to the bed during transport and unexpected power outages.

Dimensions		
Patient Surface	36" X 76"	(91cm X 193cm)
Overall Length	42" X 83"	(107cm X 211cm)
Bed Frame Clearance	6"	(15cm)
Bed Height		
Low Position	17"	(43cm)
High Position	28"	(71cm)
Bed Angles		
Max Head Elevation	75 °	
Knee Gatch	40°	
Weight		
Overall Net Weight of Bed	265lbs.	(120kg)
Max Weight Limit	500lbs.	(227kg)
Electrical		
Battery back-up(High low cycles	25 - 24V DC	
Voltage(please specify)	110V @ 60Hz or 220V @ 50Hz	
Caster & Braking System		
Caster Size	5"	(13cm)

