

# S-H-300 Hydraulic Patient Transfer Stretcher



## Key Features:

- Centrally controlled brake system
- Collapsible side rails
- Dual side foot pedals
- Retractable 5th wheel steering system
- X-Ray compatible backrest
- Shock Absorbing, non-marking bumpers
- Foley bag holder
- 3" Mattress included

## Applications:

- Emergency Room
- Recovery
- Pre-Op
- Post-Op

The Amico Care Series Stretcher collection provides caregivers the control they need in the emergency environment. The stretchers can be adjusted to various positions – raising, lowering, Trendelenburg, and reverse Trendelenburg – by actuating the hydraulic pistons. The Back Section is adjustable through the use of gas springs and the Knee Gatch can be raised and lowered with an optional manual crank.

The stretcher is made of carbon steel and painted by an electrostatic spray (powder-coat). The Base Cover is made from injection-molded ABS plastic and the locking side rails are made from a combination of steel, aluminum extrusion and plastic parts for increased durability and to facilitate easy cleaning. The Back Section includes an X-ray-compatible board.

The stretcher comes complete with an IV pole, bumper wheels and actuator dust covers to protect the hydraulic actuators.

## Length & Width

Inches / cm	83" X 33" (210.5 cm X 84 cm)
Weight Capacity	550 lb (250 kg)

## Patient Surface

Inches / cm	75.7"x 24.5"	192.4 x 62.2 cm
-------------	--------------	-----------------

## Height Range

High	33.5" (85 cm)
Low	21" (53.4 cm)

## Side Rails & Casters

Caster Diameter	8" (20.3 cm)
5th Wheel Diameter	4" (10.1 cm)

## Litter Positioning

Head	84°
Knee	30°
Trend	16.5°
Reverse Trend	16.8°

## Accessories

Mattress	Included
Removable IV Pole	Included

Amico reserves the right to change product specifications without notice.

**Warranty: Three years on parts**

Amico Corporation 85 Fulton Way, Richmond Hill, ON L4B 2N4, Canada | 71 East Industry Court, Deer Park, NY 11729, USA  
Toll Free Tel: 1.877.462.6426 | Toll Free Fax: 1.866.440.4986 Tel: 905.764.0800 | Fax: 905.764.0862