

Overbed Tables

Various specialty top styles available

Spill guard top holds up to 36 oz

Vanity drawer with mirror

Integrated accessory rail (urinal hook also available)

Shock-absorbent gas spring cylinder for easy height adjustment

Seamless cup holders

Easy, one-hand activated secondary surface

OT-F-TYYS-XU

Various base styles available

U Base H Base C Base

Lowest base on the market: 1.75"
Exposed steel; no plastic cover allows for easy cleaning

Configuration Examples


Single Top	Dual Top	Specialty	
OT-S-TYYV-XU	OT-F-TYYV-XU	OT-S-TYYN-XUJ	OT-ST-TYYN-XU-C
OT-S-TYYN-XC	OT-F-TYYS-XU	OT-S-LWBB-XU	OT-F-TXXN-XH

Customize Your Overbed Tables


Table Top Options



Accessories



Flip Table Tops*




* Flips are not available on laminate, storage, jellybean and c-base tables

Standard Top Finishes**



Standard Base Finishes



** Custom finishes are available

Storage Table Options



Storage Table Configuration Examples



Standard Dimensions

Single Table Top: 32" L x 17" W (81.3 cm x 43.2 cm)

Flip Table Top: 13.5" L x 17" W (34.3 cm x 43.2 cm)

Storage Table Top: 41.5" L x 19" W (105.4 cm x 48.3 cm)

Height Range: 28" Low - 44" High (71 cm - 111.8 cm)

U, H or C-Base: 1.75" H (4.5 cm)

Custom Dimensions Available