


Integris Gas Block Assemblies



CMS-INT-CH-O

General Specifications:

- Made To Match OEM Hill-Rom/Wittrock Integris Gas Blocks
- Includes 44" Conductive Hose with DISS Female End. Custom Lengths Available Upon Request.
- Available in Chemetron, DISS, Ohmeda, Puritan-Bennett and Oxequip Fittings.



CMS-INT-DI-O



CMS-INT-PB-A



CMS-INT-CH-V



CMS-INT-OH-O

Model Numbers
Complete Integris Gas Blocks

CMS-INT-XX-X

The "XX" Defines Gas Block Connection:

Chemetron	- CH
DISS	- DI
Ohmeda	- OH
Puritan-Bennett	- PB
Oxequip	- OX

The "X" Defines Gas:

Oxygen	- O
Medical Air	- A
MedVac	- V

CERTIFIED MEDICAL SALES